

Accelerated Surgery in 2 Steps


Available ET III/IV SS III US III




122 Taper KIT

- Surgical Kit optimized for tapered fixtures
- Simplified procedure to reduce chair time
- Optimized Drill design renders cortical drills obsolete


Soft Bone – ‘1’ Initial Drill +1



Normal Bone – ‘2’ Initial Drill +2



Hard Bone – ‘2’ Initial Drill +2



Superior Cutting Power


Optimized cutting power with improved drill design utilizing state of the art machining technology

Smooth and efficient osteotomies by reducing sequence



Bone Heating Prevention


Superior cutting power maintains low cutting temperature, despite aggressive expansion of osteotomy



ETIII Ø5.0 x 13mm, hard bone, final drill cutting temperature

Prevention

Multi-step tapered design includes pilot function – increased versatility of each drill



Smiles that last a lifetime